

GLA-1000B S

Cathode Current Cutoff Modification by Skipp May wv6f
nospam4me@juno.com



Original Dentron tubes D-50a
Replace with NOS 6LQ6/6JE6
or retrofit with EL-509 by
Svetlana

Replace RFC2-RFC5 anode
suppressors with 4 turns
Nichrome wire 1/4" dia
paralleled with 100 ohm
MF resistor, 3w approx
value 40nH

RF IN



add diode across K1



Original relay unused contacts (center) used
to short Rc (75k lwt) in transmit mode.

24 volt zener, stud mount.

